

MOPTIM<sup>®</sup>



Handheld Auto Refractometer

**easyRef Pro**

# EASY REFRACTION TEST ANYTIME, ANYWHERE

- **Shack-Hartmann Technology**

This state-of-the-art optical technology delivers unparalleled precision in eye measurement, setting a new standard for accuracy.

- **Internal Fogging System**

Designed to effectively minimize the accommodation of the eye, this feature guarantees readings of exceptional reliability and accuracy.

- **Easy Operation**

Easy to use and intuitive. Get results in less than 2 seconds per eye with a simple step: align, press, and go.



## APPLICATIONS



## SPECIFICATIONS

<b>Principle</b>	Shack-Hartmann wavefront
<b>Spherical Range</b>	-20D, +20 D
<b>Cylinder Range</b>	-12D, +12 D
<b>Axis Range</b>	1° ~ 180°
<b>Touch Screen</b>	2.80 inch
<b>Data Output</b>	Bluetooth printer

\*Specifications are subject to change due to product improvement.



Manufacturer: Shenzhen Certainn Technology Co., Ltd.  
Add: L302, Bldg.2, Skyworth Innovation Valley, No.8 Tangtou 1st Road, Tangtou Community, Shiyan Subdistrict, Bao'an District, 518108 Shenzhen, P.R.China  
Website: [www.moptim.com](http://www.moptim.com) Email: [sales@moptim.cn](mailto:sales@moptim.cn)