



# WILD FIELD FUNDUS IMAGING SYSTEM

SK-8000 SERIES



## LARGE VIEWING ANGLE

133-degree large field of view, more fundus information can be obtained in one shot  
Provide clear fundus images to help doctors accurately diagnose various fundus diseases



## MULTIPLE USES

Can be directly used for animal fundus disease inspection



## PORTABLE SIZE

Small size, easy to carry and use



# PARAMETER

	SK-8000A	SK-8000B	SK-8000C
<b>Viewing Angle</b>		133°	
<b>Minimum Pupil</b>		6mm	
<b>Collector</b>	Medical grade professional CMOS		
<b>Resolution</b>		≥10MP	
<b>Brightness Adjustment</b>		Stepless dimming	
<b>Light Source System</b>	Built-in integrated optical fiber		
<b>Shooting Method</b>	Pedal shooting, Software shooting		
<b>Shooting Mode</b>	Color photography	Color photography + FFA	
<b>Focus Method</b>	Manual focus	Auto focus	
<b>Lesion Marker</b>	Manual labeling		Automatic labeling
<b>Interface</b>	USB, DICOM		



**ChongQing Sunkingdom Medical Instrument Co., Ltd.**

Head office: Room 12, Floor 10, Block A of China Resource Center, Jiulongpo District, Chongqing 400050, China

FAX: +86 23 68102793

E-mail: [sales@cqsunkingdom.com](mailto:sales@cqsunkingdom.com)

TEL: +86 23 68643990

URL: [en.cqsunkingdom.com](http://en.cqsunkingdom.com)